

DEFINITIVE

A New Generation of  
Silicone Hydrogel Technology

DEFINITIVE

**Definitive: Leading the way**

- For those patients who need it the most, whenever oxygen transmission is important, there is only one choice - Definitive.
- No longer forced into traditional hydrogels, Definitive provides even the most challenging patients the ability to wear an appropriately healthy lens option.
- No transient surface treatment is required. Lenses made with Contamac's Definitive material are inherently wettable for outstanding on-eye comfort.
- Now your most demanding patients can be prescribed with a healthy lens choice.

CONTAMAC 

**Contamac Ltd**

Orchard House  
Bearwalden Business Park  
Wendens Ambo  
Saffron Walden  
Essex CB11 4JX UK

Tel: +44 (0)1799 544400  
Fax: +44 (0)1799 543000  
E-mail: sales@contamac.co.uk

CONTAMAC US 

**Contamac US**

806 Kimball Avenue  
Grand Junction, CO 81501

Tel: 1-(970) 242 3669  
Fax: 1-(970) 243 5501  
Toll Free: 1-866 872 6682  
E-mail: salesusa@contamac.com

contamac.com

**Definitive + Comfort + Oxygen + Design  
= Patient Satisfaction**

A New Generation of  
Silicone Hydrogel Technology



# DEFINITIVE

A New Generation of  
Silicone Hydrogel Technology




Definitive is a new Silicone Hydrogel material designed for the speciality lens market.

This unique material provides the benefits of Silicone Hydrogel Technology for all prescriptions. Now even more challenging patients have the option of a healthier, comfortable lens.

### Material Benefits...

- Unique Wettable Surface = Comfort
- Low Modulus = Comfort
- 3 x More Oxygen than Mid-Water, "Old Technology" Hydrogels = Healthy Corneas

Lens Thickness Versus On-Eye Dk/t of Definitive




Definitive + Comfort...

### Material Properties

Dk	60 Barrers
Classification	filcon II 3
Water content	74%
Modulus	≤0.39MPa
Surface	Non-treated

Higher Oxygen Transmissibility promotes corneal health. The effect of the higher O<sub>2</sub> transmissibility provided by Contamac's Definitive Silicone Hydrogel is demonstrated below:



### Contact Lens Benefits:

- Made to order, custom prescribed lenses in a Silicone Hydrogel material
- Available in a full range of designs, spherical, toric, multifocal, Keratoconus, etc.
- Excellent visual acuity
- High levels of oxygen permeability
- Highly wettable surface
- Excellent comfort

Definitive + Comfort + Oxygen...