


More opportunities to participate in the largest and fastest growing market segment.



* Sources available on file

Benefit from the choices available within the Definitive range and satisfy your practitioner's demand for Silicone Hydrogel options.

The world's leading brand now offering more choice and new opportunities

Contamac[®]

Contamac Ltd
Carlton House, Shire Hill, Saffron Walden
Essex, CB11 3AU, United Kingdom
T. +44 (0)1799 514800 F. +44 (0)1799 514801
E. sales@contamac.co.uk www.contamac.com



Contamac[®]
EMPOWERING YOUR VISION

Enhancing corneal health with a complete range of Silicone Hydrogel materials.

Whatever your requirement there is a Definitive product to suit.



Definitive range

	Water content	Oxygen permeability	Elongation to break	Modulus	Surface treatment	
<p>50</p> <p>This specialist polymer offers the opportunity for thin and durable lens designs for those labs using surface treatment techniques.</p>	50	50	290	0.6	Recommended	<p>Durable</p> <p>Thin lens designs</p> <p>Surface treatment</p>
<p>65</p> <p>Designed to retain all of the health benefits of Definitive 74 but engineered to be more durable and stable for multipack use and resistant to H₂O₂.</p>	65	62	250	0.75	Recommended	<p>Multipacks</p> <p>Parameter integrity</p> <p>Peroxide resistant</p>
<p>74</p> <p>The world's leading latheable Silicone Hydrogel: The only FDA approved material in its class.</p>	74	60	180	0.35	Not required	<p>Comfort</p> <p>Health</p> <p>Wettability</p>