

ERGO TOUCH

For a controlled delivery



“The ergonomic solution for an improved control”

Injector body design for better handling and made of crystal clear material with large fenestrations for improved visibility. Robust plunger designed with an advanced material for perfect alignment and a stopper for safe stop/start. Soft touch and anti-slip texture for safe and comfortable handling.

The ErgoTouch injector cartridges are available for a range of incision sizes from 1.8 mm to 2.6 mm.

Contamac[®]
EMPOWERING YOUR VISION

TRUSTED⁺ QUALITY
IOL MATERIALS

INNOVATION + QUALITY + CONSISTENCY

Contamac[®]

Contamac Ltd

Carlton House, Shire Hill, Saffron Walden
Essex, CB11 3AU, United Kingdom

T. +44 (0)1799 514800 F. +44 (0)1799 514801

E. sales@contamac.co.uk www.contamac.com



THE QUEEN'S AWARDS
FOR ENTERPRISE:
INTERNATIONAL TRADE
2012



THE QUEEN'S AWARDS
FOR ENTERPRISE:
INNOVATION
2016

CI26

Clear + Yellow



Hydrophilic

“A well established safe material suitable for a wide range of lens designs”

- ✚ Optimised for 2.2 mm incision “in the bag” and 1.8 mm “wound assisted”
- ≍ Water content 26 %
- ✚ Refractive index 1.460 (@ 20°C) and 1.457 (@ 35°C)
- 🔗 Tensile strength 2.5 Mpa
- ↪ Modulus 3.0 MPa

MICS22

Clear



Hydrophilic

“The next generation hydrophilic acrylic material customised for MICS IOL’s”

- ✚ Optimised for 1.8 mm incision “in the bag” and 1.4 mm incision “wound assisted”
- ≍ Water content 22 %
- ✚ Refractive index 1.465 (@ 20°C) and 1.463 (@ 35°C)
- 🔗 Tensile strength 3.9 MPa
- ↪ Modulus: 4.0 MPa
- ☀️ 10 % UV cut-off 375 nm

CI18

Clear



Hydrophilic

“Specific material properties for less conventional lens designs”

- ≍ Water content 18 %
- ✚ Refractive index 1.47 (@ 20°C)
- 🔗 Increased tensile strength
- ↪ High modulus
- ✚ Light transmission >95 % (400 – 440 nm)

PMMA

Clear



“Premium PMMA product for predictable outcome and maximum performance”

- ✚ Refractive index 1.49 (@ 20°C)