

## TECHNICAL DATA

### MICS22

#### Special Attributes

Nature: Copolymer of hydroxyethyl methacrylate and methyl methacrylate

MICS22 represents the next generation of hydrophilic acrylic material customised for micro insertion IOLs. It offers increased tensile strength, enhanced modulus, and controlled unfolding due to excellent shape memory.

MICS22's stiffness allows laboratories to design thinner IOL optics while maintaining compression forces. It's optimised for 1.8mm incision in the bag and 1.4mm incision 'wound assisted'.

#### Material Characteristics

PROPERTY	MICS22
Swell Factor at 20°C	1.11
Water Content at 20°C	22 %
Refractive Index at 20°C Hydrated	1.465
Refractive Index at 35°C Hydrated	1.463
Tensile Strength	3.9 MPa

Please note: These are typical values. Always check the certificate of compliance supplied with the goods for the actual values of the batch.

The values quoted on this data sheet are indicative of the typical values that can be obtained for the material. For the complete specification, please contact Contamac at +44 (0)1799 514800 or [sales@contamac.co.uk](mailto:sales@contamac.co.uk). Unauthorised copy or reproduction of the above information is strictly prohibited without the prior consent of Contamac Ltd.

Version: 1 - 210408

Contamac Ltd  
Carlton House, Shire Hill  
Saffron Walden, Essex  
CB11 3AU

T. +44 (0)1799 514800  
E. [sales@contamac.co.uk](mailto:sales@contamac.co.uk)  
[www.contamac.com](http://www.contamac.com)