

Accessories Brochure

For Contact Lens
Manufacturing & Handling

ContamacTM



Manufacturing 04

Laboratory Lens Handling 09

Packaging 10

Patient Lens Handling 14

Technical Information 16

To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 [0] 1799 514800

Contamac

How to read this brochure

Recommended for:



Optimum
F2
Hybrid FS



GM Advance
Contaflex
Contaflex GM3



Definitive
Unisil

[MOQ] - Minimum Order Quantity

To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 [0] 1799 514800

Contamac

Manufacturing

04



High Melt Black Wax

High melt wax for blocking Hydrophilic materials. Melting point 66°C.

Colour: Black

Quantity: 1 Tub (125 g) / 10 Sticks (20 g each) ^(MOQ)

Code: 4110047 (Tub)
4110048 (Sticks)



Low Melt Black Wax

Low melt wax for blocking RGP materials. Melting point 44-51°C.

Colour: Black

Quantity: 1 Tub (125 g) / 10 Sticks (20 g each) ^(MOQ)

Code: 4110049 (Tub)
4110050 (Sticks)



Blocking Kemdent Sticky Wax

High melt wax for blocking Hydrophilic materials. Melting point 55°C.

Colour: Yellow

Quantity: 1 Box (500 g In sticks) ^(MOQ)

Code: 4110053



To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 [0] 1799 514800

Contamac

Manufacturing



Contapol 1 Polish

High quality polishing compound suitable for hydrophilic materials, particularly those based on HEMA. It contains a uniformly dispersed alumina abrasive, suspended in an oil base, to produce excellent surfaces.

Quantity: 1 Pot (250 g) ^(MOQ)

Code: 4230001



Contapol 2 Polish

Polish specially formulated for use in high water content materials. It contains a finely graded abrasive suspended in a fine oil base, imparting excellent surface finishes.

Quantity: 1 Pot (250 g) ^(MOQ)

Code: 4230002



Contapol 3 Polish

Aqueous based premium polishing compound. Contains finely graded abrasive. Contapol 3 has a unique formula that eliminates solvent attack and enhances the wettability of the polished surface. Suitable for all RGP materials. Particularly suited to low wetting angle materials.

Quantity: 1 Bottle (1 kg) ^(MOQ)
1 Bottle (5 kg)

Code: 4230004 (1 kg)
4230005 (5 kg)



To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 (0) 1799 514800

Contamac

Manufacturing



Polishing Powder

Finely graded polishing powder for customers who wish to produce their own polishing compounds.

Quantity: 1 Tub (500 g) ^(MOQ) / 1 Drum* (10 kg)

Code: 4230006 (Tub)
4230007 (Drum)

*Available on request



Delrin Polishing Cup

High quality plastic polishing cup designed for use with our Delrin Polishing Sponge for front surface polishing.

Colour: White

Quantity: 10 Units ^(MOQ)

Code: 4110015



To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 (0) 1799 514800

Contamac

Manufacturing



Delrin Polishing Sponge

High quality foam sponge pads designed for use with our Delrin Polishing Cup.

Quantity: 100 Units ^(MOQ)
 Size: Ø 25mm x Depth 15mm
 Code: 4110034



Delrin Drum Polishing Sponge

High quality foam sponge pad designed for use in a Delrin drum polishing tool for the modification or finishing of rigid contact lenses and prosthetic lenses.

Quantity: 100 Units ^(MOQ)
 Size: Ø 7cm x Depth 7mm
 Code: 4110019



Grey Roller Sponges

Designed for use with a variety of edge polishing machines.

Colour: Grey
 Quantity: 20 Units ^(MOQ)
 Size: Ø 5.5cm x Depth 4.5cm
 (Inner circle 1cm)
 Code: 4110058



To request the Accessories
 Price List, please contact:
sales@contamac.co.uk
 +44 (0) 1799 514800

Contamac

Manufacturing

08



Microfibre Polishing Cloth

The unique construction of the micro thin fibres is particularly suited to gas permeable lens polishing. Can be used on conventional and bladder polishing systems.

Colour: Beige
Quantity: 1 Unit (61cm²) ^(MOQ)
Code: 4110005



Brushed Cotton Polishing Cloth

Universal polishing cloth, widely used throughout the industry, combines excellent results with low cost.

Colour: White
Quantity: 1 Unit (1m²) ^(MOQ)
Code: 4110006



To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 (0) 1799 514800

Contamac

Laboratory Lens Handling



Tweezers

Large lens tweezers with replaceable super soft silicone tips for increased protection when handling soft or silicone hydrogel contact lenses.

Colour: White

Quantity: 250 Units ^(MOQ)

Code: 4110059



Tweezer Tips

Replacement super soft silicone tips for white tweezers.

Quantity: 1000 Units ^(MOQ)

Code: 4110043



Laboratory Lens Handler

High quality two part green polishing suction holder designed for handling gas permeable lenses during edge polishing. Two profiles included to allow a variety of RGP designs to be held.

Colour: Green

Quantity: 2 Units ^(MOQ)

Size: Convex 8mm, Concave 9.5mm

Code: 4110039



To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 [0] 1799 514800

Contamac

Packaging



Glass Vial

Glass vial manufactured from type 1 neutral glass and fully autoclavable. Compatible with Tear Off Aluminium Seal and Siloprene Bungs. Available in two tray sizes.

Quantity: 186 Units/208 Units ^(MOQ)


Size: 5ml, 20mm

Code: 4110044 (186 Units)
4110045 (208 Units)



For technical info [186] [click here](#)

For technical info [208] [click here](#)



Siloprene Bung

Made from the finest pure siloprene, these stoppers are fully autoclavable and designed for use with our Glass Vials and both types of Crimp Top Seals. Packed without silicone lubricating oil.

Quantity: 1500 Units ^(MOQ)

Size: Ø 20mm

Code: 4110004



For technical info [click here](#)

To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 [0] 1799 514800

Contamac

Packaging

Tear Off Aluminium Seal



All aluminium seal with tear out tab that allows the complete removal of the closure prior to or after use. For use with Steel Hand Crimper Gun, Glass Vial and Siloprene Bung.

Colour: Blue, Gold, Green, Red, Silver

Quantity: 1000 Units ^(MOQ)

Size: 20mm

Code: Blue - 4110008 Red - 4110013
 Gold - 4110011 Silver - 4110014
 Green - 4110012



For technical info [click here](#)

Flip Tear Up Seal



Seal with plastic button attached to help facilitate the activation and complete removal of the seal from the container. For use with Steel Flip Off Hand Crimper Gun, Glass Vial and Siloprene Bung.

Colour: White, Blue

Quantity: 1000 Units ^(MOQ)

Size: 20mm

Code: Blue - 4110009
 White - 4110010



For technical info [click here](#)

To request the Accessories
 Price List, please contact:
sales@contamac.co.uk
 +44 [0] 1799 514800

Contamac

Packaging



Steel Hand Crimper Gun

Hand Crimper for manual crimping of Tear Off Aluminium Seals to Glass Vials.

Quantity: 1 Unit ^(MOQ)

Size: 20mm

Code: 4110031



Steel Flip Off Hand Crimper Gun

Hand Crimper for manual crimping of Aluminium Flip Tear Up Seals to Glass Vials.

Quantity: 1 Unit ^(MOQ)

Size: 20mm

Code: 4110018



To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 [0] 1799 514800

Contamac

Packaging



Double Opticase Plastic RGP Vial

Small compact RGP lens storage case features a plastic claw and two coloured lids with embossed L and R logo.

Colour: Clear with Blue and White Lids

Quantity: 250 Units ^(MOQ)

Code: 4110029

Item cannot be supplied in the EU



Double Soft Lens Mailers

High quality double mailer in clear plastic, affording excellent protection for soft lenses. This mailer can also be used for larger RGP lenses up to approximately 16mm in diameter.

Quantity: 100 Units ^(MOQ)

Code: 4110027

FDA Approved

For technical info [click here](#)



Double GP Mailers


High quality double mailers are available in clear plastic and designed for use with RGP lenses. Can be cut for use as single units.

Quantity: 250 Units ^(MOQ)

Code: 4110026

FDA Approved

For technical info [click here](#)



To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 (0) 1799 514800

Contamac

Patient Lens Handling



DMV RGP Lens Remover Ultra

High quality solid rubber lens remover designed for RGP lenses.
Supplied in individual storage cases.

Colour: Assorted

Quantity: 3 Units ^(MOQ)

Code: 4110041



Item cannot be supplied to EU



DMV Soft Lens Handler/Remover

A durable rubber removal device specifically designed for soft and silicone hydrogel lenses. Ideal for patients with limited handling ability or for those who have difficulty touching their eyes. Supplied in individual storage cases.

Quantity: 3 Units ^(MOQ)

Code: 4110033



Item cannot be supplied to EU



DMV Scleral Suction Holder

Lens inserter and remover made from high purity rubber. Designed specifically for use with scleral and other large diameter lenses.

Colour: Blue

Quantity: 3 Units ^(MOQ)

Code: 4110040



Item cannot be supplied to EU

To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 (0) 1799 514800

Contamac

Patient Lens Handling



RGP Lens Inserter/Remover

A lower cost alternative RGP lens inserter and remover, hollow for increased suction. Supplied in individual storage cases.

Colour: Green

Quantity: 10 units ^(MOQ)

Code: 4110037

Item cannot be supplied in the EU



Soft Lens Remover

A lower cost alternative solid soft lens remover. Supplied in individual storage cases.

Colour: Orange

Quantity: 10 Units ^(MOQ)

Code: 4110038

Item cannot be supplied in the EU




To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 [0] 1799 514800

Contamac

Technical Information


16

Glass Vial (186) 5ml



d1	22,25	+0,25	-0,25	h3	7,75	+0,5	-0,5	s1	1,1	+0,05	-0,05
d2	19,9	+0,2	-0,2	h21	3,0	+0,2	-0,2	r21	0,5	+0,2	-0,2
d3	16,05	max.		c21*	30,0	+5,0	-5,0	r23	0,5	ca.	
d4	12,75	+0,15	-0,15	c22*	5,0	+2,0	-2,0	t	0,4	+0,3	-0,3
h1	38,0	+0,5	-0,5	c23*	45,0	ca.		a	0,66		
h2	24,5	min.		c24*	45,0	ca.					

Siloprene Bung 20mm






MATERIAL: Type 1, Clear Neutral Glass

Richtmasse Dimensions Indicatives Standard Dimnsions	19.3	13.1	3.5	3.9	8.4
	øA±0.1	øB±0.1	D±0.3	E±0.1	H±0.4

Technical Information




18

Tear Off Aluminium Seal



[Product Info](#)

Flip Tear Up Seal




To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 [0] 1799 514800

Contamac

Technical Information

19

Double GP Mailer



[Product Info](#)


To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 [0] 1799 514800

Contamac

Technical Information

20

Double Soft Lens Mailers



To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 [0] 1799 514800

To request the Accessories
Price List, please contact:
sales@contamac.co.uk
+44 (0) 1799 514800

Contamac[™]

